



# Drug of Abuse (DOA)

## One-Step Rapid Test



Accuracy  
Fast  
Excellent

Accuracy more than 99%  
Easy to read, Easy to use  
Up to 14 drugs test at once  
Format strip, cassette,  
multi-strip and cup  
False warranty

# Drug of Abuse Test

Specimen

✓ Urine

Description	Sensitivity (ng/ml)
Amphetamine	1000
Barbiture	300
Benzodiazepine	300
Cocaine	300
Marijuana (THC)	50
MDMA / Ecstasy	500
Methamphetamine	1000
Morphine (Opiate)	300
DOA-3 Panel Test (AMP/MOR/THC)	1000/300/50
DOA-5 Panel Test (AMP/MOR/THC/MET/BZD)	1000/300/50/ 1000/300
DOA-6 Panel Test (AMP/MOR/THC/MET/BZD/COC)	1000/300/50/ 1000/300/300
DOA Panel Test Cup up to 14 drugs, split type	
DOA Panel Test Cup up to 14 drugs, on-site type	

## Product information

**Accuracy** : 99,86%, **Specificity** : 99,93%, **Sensitivity** : 99,80%

### Production (Shelf Life) :

24 months

### Product Formats

Strip, Cassette, Multi Strip, Cup

### Packaging

50 test strips per box  
25 test cassettes per box  
25 test multi strips per box  
1 cup per pcs

## Specification

### Method

Immuno-chromatography assay for qualitative detection

### Reading Time

5 minutes

### Specimen

Urine

### Size

Strip : 4mm x 75mm  
Device : 25mm x 70mm x 5mm  
Multi 3 : 31mm x 80mm x 5mm  
Multi 5 : 43mm x 80mm x 5mm  
Multi 6 : 50mm x 80mm x 5mm

ISO 13485

GMP CE