






## RT - 1904CV for veterinary

**Semi Auto Chemistry Analyzer** 



**Menu Selection** 

## **System Spesification**

Principle Absorbance

Photometric range -0.500 ~ 3.000 Abs

Resolution 0.001Abs (displayed), 0.001Abs (Calculated)

Light Source Halogen Lamp

Wavelengths 340, 405, 500, 546, 578, 620, 670nm,1 more filter optional

Wavelength Accuracy ± 1 nm

Band Width 8 nm

Aspiration Volume 100 ~ 3000 ul

Carryover <0.1 %

Memory 54 test parameters, up to 500 sample results

Interface RS-232, USB, SD card interface
Flow Cell 25ul Metal-quartz flow cell

Cuvette Mode 12.5mm x 12.5mm cuvette (optional)

Temperature Control 25°C, 30°C, 37°C, ±0.1°C and ambient temperature

Display 7" color LCD (640x240, 256 colors)
Input USB mouse and external keyboard

Output Build-in printer or external printer (optional)

Net Weight 8 kgs

Dimension L x W x H (mm)  $460 \times 330 \times 140$ 

Power Supply AC 110-220V ±10%, 50-60Hz



Easy Windows operation system with mouse.

Large LCD Color display.

Support up to 12 kinds of animals: 8 kinds fixed, 4 kinds open.

0000

A ayto

34 most common tests fixed, 20 more can be defined by user.

Reagent open system supports flow cell and cuvette mode,

Real time curve monitoring enhances the reliability of result.

Powerful QC Function : Westguard multi-rule, Levey Jennings plot, 2 controls per test.

Auto self-check when power on.

With lamp saving and wake-up function.

Easy software upgrade by SD card, support USB mouse and keyboard.

Compact design enables USER-DIY maintenance.

More than 500 samples results data can be stored.



Rayto Life and Analytical Sciences Co., Ltd.

2nd Floor, 5th Building, Software Park, NO.2,Gaoxinzhong 3rd Road,High-tech Industrial Park, Nanshan, Shenzhen 518057,P.R.China Tel:+86(755)86168996 Fax:+86(755)8616879

Email: info@rayto.com Http://www.rayto.com